


## Type Ca ECG Monitor Cables




### Ca-type ECG monitor main cable

Product No: 0202230-I  
Model No: CaD-12Pu/IEC5/Q



### Ca-type ECG monitor lead wire



### Electrodes for Ca-type ECG monitor lead wire



The electrodes for Ca-type ECG monitor lead wire are open-and-close electrodes.

The shape of the electrode net tail is designed to be curved, therefore the wire will not be bent in use. The service life and the quality of the signal are better guaranteed.

Patent No. : 201620979821.9

### Application demonstration

