





ECG/EKG Cables & Lead wires

Compatible Brand	Product Picture and NO.	Product Description	Application	Drawing
NIHON KOHDEN	0201002 	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ3.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201002-I Product no for AHA: 0201002-A Spec&Model: NK-B-15Pa-B	Compatible with NIHON KOHDEN ECG-6511, 6151, 6551, 8110P/K 6353, 6403, 6504 Zhonglianhaiuo XD-106,XD-104B/D ect ECG(Universal)and Dongjiang ECG	
NIHON KOHDEN	0201003  Patent NO.: ZL200930331335.1	Plug for Main Cable:Special 15PIN plug Needle dimension: Φ3.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201003-I Product no for AHA: 0201003-A Spec&Model: BD-15Pb-B	Compatible with ECG9000,9010,1250,/9020, 9020K, 9130, 9100,9103K,9130P,9552P, 9020P,9022,9110,9620,9620P,9132, 92C 8420 series ECG (except 9320 9522P,1550P)	
FUKUDA	0201037 	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ3.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 950mm Length of Limbs: 1650mm & 1850mm Product no for IEC: 0201037-I Product no for AHA: 0201037-A Spec&Model: GDF-15PY-B	Compatible with Fukuda FX-2100,2150,4010, FX2111,2155,3010,7102,7202, 7402,FCP-7411 series ECG	
PHILIPS/HP	0201008  Patent NO.: ZL200930331335.1	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ3.0mm Plug Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201008-I Product no for AHA: 0201008-A Spec&Model: BaD-15PH-M	Compatible with Philips/HP M3703C,M2462A, M1770A PAGEWriter 200i, M1771A PAGEWriter 200, M1772A PAGEWriter 100	
GE (MARQUETTE)	0201014 	Plug for Main Cable: 十组合病人导联线 Needle dimension: Φ4.0mm banana Plug Length of Chest Lead Cent Cable: 800mm Length of Limbs: 1150 Product no for IEC: 0201014-I Product no for AHA: 0201014-A	Compatible with GE (Marquette) Mac-500、 Mac-800、 Mac-1200, Mac-5000 ect ECG	
GE (MARQUETTE)	0201027 & 0201014  Patent NO.: ZL 200720027004.4 ZL 200520086392.4	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ4.0mm banana Plug Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 800mm Length of Limbs: 1150 Product no for IEC: 0201027-I & 0201014-I Product no for AHA: 0201027-A & 0201014-A Spec & Model: KD-15Pid/KLL-M	Compatible with GE (Marquette) Mac-500、 Mac-800、 Mac-1200, 5000 ect ECG	
GE (MARQUETTE)	0201085  Patent NO.: ZL200810249520.0	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ4.0mm banana Plug Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm&1250mm Product no for IEC: 0201085-I Product no for AHA: 0201085-A Spec & Model: TD-15Pid-J/AHA10/Q	Compatible with GE (Marquette) Mac-500、 Mac-800、 Mac-1200, Mac-5000 ect ECG Special recommendation	
GE (MARQUETTE)	0201034 & 0201014  Patent NO.: ZL 200520086392.4	Plug for Main Cable: Special 11PIN plug Needle dimension:Φ 4.0mm banana Plug Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 800mm Length of Limbs: 1150mm Product no for IEC: 0201034-I & 0201014-I Product no for AHA: 0201034-A & 0201014-A Spec & Model: KD-11Pd/KLL-M	GE(Marquette)Mac Eagle,Solar, Dash Monitors,Tram, MAC-Lab Cath Lab System ect ECG	
GE (MARQUETTE)	0201027 & 0201015  Patent NO.: ZL 200720027004.4 ZL 200520086392.4	Plug for Main Cable: Special 15PIN plug Electrode Part: Button Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 800mm Length of Limbs: 1150mm Product no for IEC: 0201027-I & 0201015-I Product no for AHA: 0201027-A & 0201015-A Spec&Model: KD-15Pid / KLL-N	Compatible with GE (Marquette) dynamic treadmill Mac-case ECG	
SCHILLER	02011156 	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ 4.0mm banana Plug Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 770mm Length of Limbs: 1250mm Product no for IEC: 0201038-I Product no for AHA: 0201038-A Spec & Model: MDR-15Ph/9-M	(Schiller) Cardiovit AT-1 、 AT-2 , AT-102 , CS-200 , AT-10 / PLUS ECG , AT-104 PC, AT-110 , AT-4 , Meilin : Eii250	
SIEMENS	0201089  Patent NO.: ZL200930331335.1	Plug for Main Cable: Special 15PIN plug Needle dimension: Φ3.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201089-I Product no for AHA: 0201037-A Spec&Model: BaD-Siemens15P-J	Compatible with Siemens ECG machines	
Compatible	0201017 	Plug for Main Cable: Special 15PIN plug Electrode Part: Button Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201017-I Product no for AHA: 0201017-A Spec & Model: BD-15Pb-N	Many kinds of domestic ECG & ECG workstation , With defend defibrillation	
Compatible	0201018  Patent NO.: ZL 200520086388.8 ZL 200810249520.0	Plug for Main Cable: Special 15PIN plug Electrode Part: Combined European Clamp Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201018-I Product no for AHA: 0201018-A Spec & Model: TD-15Pk-MI	Many kinds of domestic ECG & ECG workstation , With defend defibrillation	
Compatible	0201021 	Plug for Main Cable: Special 15PIN plug Electrode Part: Φ4.0mm banana Plug Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm & 1250mm Product no for IEC: 0201021-I Product no for AHA: 0201021-A Spec & Model: Qa-15Pk-Jm	Many kinds of domestic ECG & ECG workstation , With defend defibrillation	
Compatible	0201053  Patent NO.: ZL 200930331335.1	Plug for Main Cable: Special 15PIN plug Electrode Part: Φ 4.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm&1250mm Product no for IEC: 0201053-I Product no for AHA: 0201053-A Spec & Model: BaD-15PK-J	Many kinds of domestic ECG & ECG workstation , With defend defibrillation	
Compatible	0201084 	Plug for Main Cable: Special 15PIN plug Electrode Part: Φ 3.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm&1250mm Product no for IEC: 0201084-I Product no for AHA: 0201084-A Spec & Model: Qa-15Pk-R	Many kinds of Mini-ECG or ECG detected machine	
Compatible	0201088 	Plug for Main Cable: Special 15PIN plug Electrode Part: Button Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm&1250mm Product no for IEC: 0201088-I Product no for AHA: 0201088-A Spec & Model: Qa-15Pk-Da	Many kinds of Mini-ECG or ECG detected machine	
Compatible	0201090 	Plug for Main Cable: Special 15PIN plug Electrode Part: Φ 3.0mm Needle Length of Main Cable: 2100mm Length of Chest Lead Cent Cable: 670mm Length of Limbs: 1150mm&1250mm Product no for IEC: 0201090-I Product no for AHA: 0201090-A Spec & Model: BaD-15PK-BK/CK	Many kinds of domestic ECG & ECG workstation , With defend defibrillation	