



ECG Simulator For IOS

Product Introduction

Wireless Analog Signal Generator is a machine by handset, can be widely used on Multi-parameter Monitor, ECG, and Holter ect.

The Analog Signal Generator can simulate 12-lead ECG signals, 16 kinds of ST segment-change signals, 20 kinds of arrhythmia signals, and 18 kinds of pacemaker signals.

Product Features

Realize App records stored for various kinds of ECG software.

Transfer and receive the signals via Bluetooth.

Realize the App upgrade.

Connect with various kinds of ECG Machines, Monitor etc.

Be compatible with some ECG functions of Fluke Patient Simulator.


With IOS system, Bluetooth 3.0 Agreement



IOS
iphone & ipad
APP



Choice of different models



Product Display Effect

