
















































# Patient Monitor Cables

Compatible Brand	Product Picture and NO.	Product Description	Application	Drawing
Compatible	 Patent NO.: ZL 200720027017.1	Product Name: Three-lead Monitor Cable Spec&Model: ND-AMP-N Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Style of joint with electrodes: Button Product no for IEC: 0202003-I Product no for AHA: 0202003-A	Domestic and overseas Monitor machines With defend defibrillation	
	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec & Model: ND-Lemon 6P-N Plug: Lemon 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Style of joint with electrodes: Button Product no for IEC: 0202008-I Product no for AHA: 0202008-A	Domestic and overseas Monitor machines With defend defibrillation	
	 Patent NO.: ZL 200720027017.1	Product Name: Three-lead Monitor Cable Spec & Model: ND-AMP-N Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Style of joint with electrodes: Button Product no for IEC: 0202090-I Product no for AHA: 0202090-A	Domestic and overseas Monitor machines With defend defibrillation	
	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Three-lead Monitor Cable Spec & Model: N-6Po-MI Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Electrode Part: Combined European Clamp Product no for IEC: 0202093-I Product no for AHA: 0202093-A	Domestic and overseas Monitor machines With defend defibrillation	
	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec&Model: ND-AMP-N Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Style of joint with electrodes: Button Product no for IEC: 0202091-I Product no for AHA: 0202091-A	Domestic and overseas Monitor machines With defend defibrillation	
Compatible	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Three-lead Monitor Cable Spec & Model: N-6Po-MI Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Electrode Part: Combined European Clamp Product no for IEC: 0202094-I Product no for AHA: 0202094-A	Domestic and overseas Monitor machines With defend defibrillation	
	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec&Model: ND-AMP-N Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Style of joint with electrodes: Button Product no for IEC: 0202091-I Product no for AHA: 0202091-A	Domestic and overseas Monitor machines With defend defibrillation	
	 Patent NO.: ZL 200720027017.1	Product Name: Three-lead Monitor Cable Spec & Model: N-6Po-MI Plug: AMP 6PIN Length of main cable: 2700mm Length of cent cable: 900mm Electrode Part: Combined European Clamp Product no for IEC: 0202094-I Product no for AHA: 0202094-A	Domestic and overseas Monitor machines With defend defibrillation	
GE (Marquette)	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Three-lead Monitor Cable Spec & Model: KLL-5N Plug: MAC special interface Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Button Product no for IEC: 0202038-I Product no for AHA: 0202038-A	Compatible with Mac-Dash 2000(three-lead, five-lead), Mac-Dash3000, 4000 (five-lead), at the same time compatible with Eagle 4000 series, Solar 8000 series etc. Monitor machines.	
	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Five-lead Monitor Cent Cable Spec & Model: KLL-MI Plug: MAC special interface Length of cent cable: 950mm Style of joint with electrodes: Button Product no for IEC: 0202036-I Product no for AHA: 0202036-A	Compatible with Mac-Dash 2000(three-lead, five-lead), Mac-Dash200, 2500, 3000, 4000(five-lead), at the same time compatible with Eagle 1000, 3000 series, Solar 7000, 8000 series etc. Monitor machines.	
GE (Marquette)	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Five-lead Monitor Main Cable Spec & Model: KI-11P Plug: MAC 11PIN Length of main cable: 2700mm Product no for IEC: 0202064-I Product no for AHA: 0202064-A	Compatible with Mac-Dash 2000, Dash 3000, Eagle series, SOLAR series, 7000 series, Tram modular series, MAC-Lab, Cath-Lab With Resist elec-generator, defend defibrillation	
	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Five-lead Monitor Cable Spec & Model: ND-11Pa-MI Plug: MAC 11PIN Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Combined European clamp Product no for IEC: 0202109-I Product no for AHA: 0202109-A		
PHILIPS/HP	 Patent NO.: ZL 200720027017.1	Product Name: Three-lead Monitor Cable Spec & Model: ND-12Pt-N Plug: HP 12 PIN Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Button Product no for IEC: 0202098-I Product no for AHA: 0202098-A	Philips/HP 78352, 7834C 78834C, M1001A/B M1002A/B M1722A/B M1723A/B M2475B M3001A M4735A 43100A 43110A 43120A 78354C 862474 C3 M1165A M1166A M1167A M1175A M1205A	
	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec&Model: ND-12Pt-N Plug: HP 12 PIN Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Button Product no for IEC: 0202099-I Product no for AHA: 0202099-A		
PHILIPS/HP	 Patent NO.: ZL 200720027017.1	Product Name: Three-lead Monitor Cable Spec & Model: ND-12Pt-MI Plug: HP 12 PIN Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Combined European clamp Product no for IEC: 0202106-I Product no for AHA: 0202106-A		
	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec & Model: ND-12Pt-MI Plug: HP 12 PIN Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Combined European clamp Product no for IEC: 0202107-I Product no for AHA: 0202107-A	M1264A M1275A Transport, M1280A, M1722A/B, M2601A, M3046A M2/M3/M4	
SPACELABS	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Main Cable Spec & Model: WD-17Pa-O Plug: Spacelabs special 17PIN Length of main cable: 2600mm Plug-pin: Special interface Product no for IEC: 0202069-I Product no for AHA: 0202069-A	Spacelabs 1050, 1600, 1700, 90369, 90496	
SPACELABS	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cent Cable Spec & Model: 2Pa-N Plug-pin: Spacelabs special interface Length of cent cable: 950mm Style of joint with electrodes: Button Product no for IEC: 0202028-I Product no for AHA: 0202028-A	Spacelabs 1050, 1600, 1700, 90369, 90496	
SPACELABS	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec & Model: N-17P-N Plug: Spacelabs special 17PIN Length of main cable: 2700mm Length of cent cable: 950mm Style of joint with electrodes: Button Product no for IEC: 0202101-I Product no for AHA: 0202101-A	Spacelabs 1050, 1600, 1700, 90369, 90496	
PRIMEDIC	 Patent NO.: ZL 200520086388.2 ZL 200520086388.8	Product Name: Four-lead Monitor Cable Spec & Model: L-miniR-Lemo8p-MID Plug: Lemo8PIN Length of main cable: 2700mm Length of cent cable: 950mm Electrode Part: Combined European Clamp Product no for IEC: 0202072-I Product no for AHA: 0202072-A	Apply to PRIMEDIC, Remove-tremble Monitor EC01, DM1/DM3, DM10/DM30, With defend defibrillation	
NIHON KOHDEN	 Patent NO.: ZL 200720027017.1 ZL 200520086388.8	Product Name: Three-lead Monitor Cable Spec & Model: ND-AMP-MI Plug: AMP 6PIN Length of main cable: 2750mm Length of cent cable: 900mm Style of joint with electrodes: Combined European clamp Product no for IEC: 0202129-I Product no for AHA: 0202129-A	Compatible with Nihon Kohden 7531 series ECG monitor machine, With defend defibrillation	
	 Patent NO.: ZL 200720027003.X ZL 200520086388.8 Patent Application NO.: 2010301342719	Product Name: ECG Monitor Cable Spec & Model: Qa-12Pu-Mi Plug: 12PIN Length of main cable: 2050mm Length of chest cent cable: 700mm Length of Limb cable: 1000mm Electrode Part: Combined European Clamp Product no for IEC: 0202130-I Product no for AHA: 0202130-A	Apply to Schiller DG5000, With defend defibrillation	
Artema	 Patent NO.: ZL 200720027017.1	Product Name: Five-lead Monitor Cable Spec & Model: 10P-N/AHA5 Plug: 10P Length of main cable: 2750mm Length of cent cable: 900mm Electrode Part: ø4.0mm banana Plug Product no for IEC: 0202135-I Product no for AHA: 0202135-A	Ohmeda, Artema	