





Sleepmate: Respiratory Inductive Belt Adult (Disposable)



1200mm



1500mm

Product Instruction

The elastic fabric weave, good air permeability. Length can be customized according to customer requirements. Independent packing is convenient to use. Signal derived by using a common electrode button, easy matching signal transmission line. Applicable to the collection of thoracic and abdominal respiratory motion waveform.

Special cables



Product Instruction

Respiratory Inductive Belt Adult using a dedicated medical signal wire, the signal transmission and the bending resistant performance. According to customer requirements supporting lead wire and machine commutator. Double cable design, cable branch junction can be freely adjusted branch mouth position.