



My Heart™ DISPOSABLE ECG ELECTRODE

REF.DCODE:0801170

BASIC PRODUCT INFORMATION

Shape:	Oval
Size :	36*42mm
Adhesive area:	873mm ²
Gel size:	314mm ²
Sensor size:	78.5mm ²
Stud size:	Ø3.9mm
Foam thickness:	1.0mm

MATERIALS INFORMATION

Substrate Material:	White Polyethylene Foam
Adhesive:	Acrylic, Medical Grade
Adhesive Strength:	Regular
Gel Type:	Conductive Solid Gel
Sensor:	Molded of 20% Glass-filled ABS coated with Ag/AgCl
Adaptor / Connector (Stud):	Stainless Steel
Liner:	Polyethylene film (PE), one side siliconized

ELECTRICAL PERFORMANCE

(ANSI / AAMI EC12)

ACZ impedance @10 Hz:	Less than 100 mu A (peak to peak)
Average AC impedance:	2000Ω
DC offset voltage:	≤100MV
Internal Noise:	≤150 μV (peak to

Simulated defibrillation recovery performance:

- Discharge in the 5th second electrode to voltage:	≤ 100MV
- Average rate of change every 10 second:	≤ 1MV/S
Simulated defibrillation Electrochemical Impedance Spectroscopy:	≤ 3KΩ

SPECIAL FEATURES

MRI Conditional:	No
X-ray Translucence:	No

BIOCOMPATIBILITY

TEST ACC. TO DIN ISO 10993:	PASSED
Latex Content:	NO

REGULATORY

CE Mark according to MDD 93/42/EEC

SHELF LIFE

Product Shelf Life:	24months (in accordance with storage guidelines)
---------------------	--

CATAGORY

Resting ECG:	Yes
Stress ECG:	Yes
Holter:	Yes (suggest change after 72 hours)

PACKAGING

Product Packaging Materia:	Foil Paper
Packing:	50pcs/Pouch
Pouch Size (L x W):	30*20cm
Transport Packaging:	Carton/5000pcs
Carton Size (L x W x H):	**cm
Weight of the Carton:	**kgs

Disclaimer: Specifications and other information presented in this document are for information purpose only and do not provide any recommendation or advise of any kind.
Data subject to changes without notice.

My heart, Qingdao Bright with logo and ™ marked brands are trademarks of Qingdao Bright Medical Manufacturing Co.,Ltd..
© 2018 Qingdao Bright Medical Manufacturing Co.,Ltd. All rights reserved.