



# Clamp Electrodes

| Product NO. | Picture   | Description  | Application  |
|-------------|---|--|--|
| 0102004     |                             | Clamp's length: 115mm<br>(red yellow green black)<br>With plug-pin: 3.0mm<br>Material: ABS、Cu / Ni<br>Apply for: Adult<br>Spec & Model: B30-115*4/IEC4   | Nihon Kohden, Fukuda, KENZ, Burdick, Dongjiang, HP etc, ECG and ECG workstation.   |
| 0102005     |                            | Clamp's length: 115mm<br>(red yellow green black)<br>With plug-pin: 3.0mm<br>Material: ABS、Cu / Ni<br>Apply for: Adult<br>Spec & Model: B30-115*4/IEC4   | Nihon Kohden, Fukuda, KENZ, Burdick, Dongjiang, HP etc, ECG and ECG workstation.   |
| 0102026     | <br>专利号: ZL 201020253936.2 | Red、Yellow clamp's length: 115mm<br>Green、Black clamp's length: 142 mm<br>Two Big、Two Small<br>With plug-pin: 3.0mm/ 4.0mm<br>Material: ABS、Cu / Ni<br>Apply for: Adult<br>Spec & Model: B30/40-115*2+142*2/IEC4 | Nihon Kohden, Fukuda, KENZ, Burdick, Dongjiang, HP etc, ECG and ECG workstation, matched with 3.0mm/ 4.0mm needle cable. |
| 0102007     |                            | Clamp's length: 140mm<br>4 pcs of blue clamps<br>With plug-pin: 4.0mm<br>Material: nylon、Cu/Agcl<br>Apply for : Adult<br>Spec & Model: A26-A30-21/G6+2   | HP/Philip, GE/Marquette, Siemens, Schiller ect ECG (IEC)   |
| 0102008     |                            | Clamp's length: 140mm<br>Colorful clamp electrodes<br>With plug-pin: 4.0mm<br>Material: nylon、Cu / Agcl<br>Apply for: Adult<br>Spec & Model: C40-137*4/IEC4  | HP/Philip, GE/Marquette, Siemens, Schiller ect ECG (IEC)   |
| 0102009     |                            | Clamp's length: 89mm<br>Colorful clamp electrodes<br>With plug-pin: 4.0mm<br>Material: nylon、Cu / Agcl<br>Apply for: Children<br>Spec & Model: D40-89*4/IEC4   | HP/Philip, GE/Marquette, Siemens, Schiller ect ECG (IEC)   |
| 0102014     |                            | Clamp's length: 140mm<br>Colorful clamp electrodes<br>With plug-pin: 3.0mm / 4.0mm<br>Material: nylon、Cu / Agcl<br>Apply for: Adult  | HP/Philip, GE/Marquette, Siemens, Schiller ect ECG (IEC)   |