




Qingdao Bright



Disposable Sterile Subdermal Electrodes




Product Instruction

Probe with medical grade stainless steel wire, with a special process from grinding, can reduce the patient puncture bleeding and pain.

Using varied different colors on different parts of coding, easy to distinguish.

Using the United States DuPont Co dialysis paper making pocket packaging, through the epoxy ethane sterilization, pocket packaging with disinfection of ethylene oxide sterilization indication, change color after Ethylene oxide sterilization



Product Instruction

Disposable Sterile Subdermal Needle Electrodes (twisted-pair) could effectively reduce the electrical interference and noise of the effects of signal. More varieties of different specification; needle length is divided into 13mm, 20mm, 30mm; the length of the wire is divided into 1000mm, 1500mm, 2000mm, 2500mm, can produce according to customer requirements.

Can be developing new products according to the special requirements of customers.



Qingdao Bright Medical Optics Co., Ltd.

Add: NO.1 TianKang Rd, Chengyang Dist, Qingdao, China

Tel: 0086-532-67767700

FAX: +86-532-8380 1122ext 1

E-mail: export@bright-med.com

URL: [Http://www.bright-med.com](http://www.bright-med.com)