



EEG/BCI/ICG Supporting Solutions

Product introduction:


Size: L/S/M

Compatible models: compatible with various EEG machine, used for transmitting bioelectricity signals

- Available in various sizes and specification, suitable for different users.
- Widely compatible with various imported and domestic EEG equipments



EEG/BCI/ICG Supporting Solutions



EEG cap matching products



Name of product: Lead wire for EEG cap
 Device connector: standard socket $\phi 1.50-\phi 1.60\text{mm}$
 Wire length: 1500mm
 Wire diameter: $\phi 1.1, \phi 1.2$
 Electrode connector: electrode contactor
 Electrode material: Ag/AgCl
 Compatible with: designated model of EEG caps only