




Sleepmate: Respiratory Inductive Belt Adult

Sleepmate:Respiratory Inductive Belt Adult(Reusable)



Product Instruction

The electrode division supporting various fixed modes;through the bandage is fixed,through medical permeable adhesive tape fixed.
 Reproducibility of thoracoabdominal sensing zone.
 The elastic fabric weave,good air permeability.
 Can freely adjust the length,adjust the range of 600-1800mm,can be customized according to customer requirements.
 Signal derived by using a common electrode button,easy matching signal transmission line.
 Applicable to the collection of thoracic and abdominal respiratory motion waveform.



Package mode(bag)

Special cables



Product Instruction

Respiratory Inductive Belt Adult using a dedicated medical signal wire,the signal transmission and the bending resistant performance.
 According to customer requirements supporting lead wire and machine cummutator.
 Double cable design,cable branch junction can be freely adjusted respiratory mouth position.